

TRADE **"PRENTIF"** MARK

THE HALL-MARK OF EXCELLENCE



"PRENTIF" COMPOUND

For use with the "Prentif" Patent Cavity Rim Pessary
(and all rubber Pessaries), Contraceptive Sponges, etc.

PER 2/6 TUBE

3 TUBES 7/-

The standard Compound is prepared with Quinine as the
Spermaticide but it is also supplied with Chinosol or Lactic
if required.

TRADE **"PRENTIF"** MARK

LUBRICATING JELLY

The ideal lubricant for condoms and sheaths. The "Prentif"
Lubricant contains no grease and will not harm or perish
rubber.

PER 9d TUBE

OBTAINABLE FROM:

W. GEORGE'S LTD.

10, GREEN ST. LEICESTER SQ. LONDON, W.C.2.

Mail orders should also be sent to above

or from branches

125a, QUEEN'S RD. S.E.15.

21, GREEN ST. LEICESTER SQ. W.C.2.

39, MONTPELIER VALE, BLACKHEATH, S.E.3.

1, BURNT ASH ROAD, LEE, S.E.12.

Full and Illustrated Price List will be sent on application -
post free and under private cover.

Trade **"PRENTIF"** Mark

WASHABLE SHEATH

When not in use, the Sheath should be kept unrolled in the box. Before use, it should be rolled up (between the fingers and thumbs) on to its rim ; the top should then be squeezed to exclude the air and the sheath unrolled into position, a space being left at the top. "PRENTIF" Lubricating Jelly should be applied to the outside of the sheath, from the top to about two inches down.

If the sheath is not being washed immediately after use, it should be placed in water. For washing, use only tepid water, dry both sides thoroughly with a soft cloth and powder the sheath inside and out with best quality French Chalk or an unmedicated talcum powder.

IF ANYTHING GREASY-FOR EXAMPLE, VASELINE
OR COCOA-BUTTER - IS ALLOWED TO COME
INTO CONTACT WITH THE SHEATH, ITS LIFE
WILL BE CONSIDERABLY SHORTENED

"PRENTIF" Jelly - The perfect lubricant. Non-greasy and cannot harm the rubber. PRICE 9d. per tube.

Sole Manufacturers - BCM/PRENTIF, London, W.C 1.

OBTAINABLE FROM ALL CHEMISTS

PATENT

PRELIMINARY

DEAL PREVIEW

EXHIBIT

CONFIDENTIAL

AND

STYLING

SYSTEM