

Introducing
the...

DURAFLEX

A Requisite in the Treatment
of Cases Presenting Patho-
logical Changes in the Birth
Canal.




(Pat. applied for)

Distributed to members of the medical
profession or at their direction

by


DUREX PRODUCTS, INC.
156 FIFTH AVE. NEW YORK, N.Y.

raf
ec
ace
sai
ctu



Cystocele condition.



One type of case in which the Duraflex will prove helpful.



Duraflex type lifting anterior vaginal wall with cystocele.


Ineffectiveness of round spring vaginal diaphragm — in case of cystocele.

(This is especially true if weak spring diaphragm is employed.)



Simple and effective introduction of the Duraflex with aid of the Mayer Pessary Director.

Duraflex Corrective and/or Fitting Rings



Size D —
(approx. 100 mm long)



Size G —
(approx. 85 mm long)

Duraflex Corrective Contraceptive Pessaries

(All actual size)




← Size D
(approx.
100 mm long)



Size E →
(approx.
95 mm long)




← Size G
(approx.
85 mm long)




Size F →
(approx.
90 mm long)

Duraflex Corrective and/or Fitting Rings



Size E —
(approx. 95 mm long)



Size F —
(approx. 90 mm long)