


THE SCHMITZ

DIAPHRAGM INTRODUCER

*AN IMPORTANT
Contribution to
the Technique of Fit-
ting Diaphragmatic
Pessaries.*




Holland-Rantos Co., Inc.
37 East 18th Street
New York, N. Y.



THIS professional instrument—compact and hygienic—is a simple yet ingenious device for placing the diaphragm in correct position without manual manipulation of the pessary.

☛ The illustrations picture clearly its ease of use. The diaphragm, symmetrically stretched on the tip of the instrument,



is deftly guided to its anchorage between the posterior cul-de-sac behind the cervix and under the pubic bone. A slight pressure of the finger on the Introducer's plunger gently releases the diaphragm, as the instrument is withdrawn. ☛ For fitting purposes the Introducer is provided with a series of tips to accommodate diaphragms ranging in size from 55 to 85 millimeters. For use by the patient only one tip, corresponding to the size of diaphragm prescribed by the physician, is needed.



☛ The Introducer which is made of Bakelite, and thus easily sterilized, is put up in a box with four tips (65-70-75-80) for the physician's use in fitting. Tips sold separately.



HOLLAND-RANTOS COMPANY, Inc., Distributor